



AdCap® is available with silica gel, molecular sieve and blends.

1.0 g

∅ 14.0 × 17.6 mm

2.0 g

∅ 19.4 × 15.8 mm

All dimensions are approximate. Other sizes on request.

Any user of this product is responsible for determining the suitability of Sanner products for each application.

109.02.07.2015

SANNER GmbH
Schillerstraße 76
64625 Bensheim
Germany
European Sales Office
Phone +49(0)6251 938-0
Fax +49(0)6251 74672
sales@sanner-group.com
www.sanner-group.com


Sanner of America, Inc.
2400 Arbor Tech Drive
Hebron, Kentucky 41048
USA
U.S. Sales Office
Phone +1 859 795 1955
Fax +1 859 795 1954
sales@sanner-usa.com
www.sanner-group.com

Sanner Pharmaceutical & Medical Packaging Materials Co., Ltd.
Asian Sales Office
No. 819 Nanjing Road (W)
Zhongchuang Building,
Room 1308 Shanghai 200041
China
Phone +86-21-52130427
Fax +86-21-52130422
Production
No. 666 Nanbang Road
Kunshan Jiangsu Province
215300 China
info@sanner.com.cn
www.sanner.com.cn

Sanner Hungária Kft
Haraszti út 38/A
1239 Budapest
Hungary
European Sales Office
Phone +49(0)6251 938-0
Fax +49(0)6251 74672
sales@sanner-group.com
www.sanner-group.com

Sanner Indonesia
Indonesia Sales Office
Talavera Office Park, 28th Floor
JL.TB. Simatupang Kav. 22-26
Jakarta 12430 Indonesia
Phone +62(21) 75999889
Fax +62(21) 75999888
info@sanner.com.cn
www.sanner.com.cn

Sanner India
c/o IndiGate
100/5 Ambegaon
Pune – 411 046
India
Phone +91(020) 39392211
Fax +91(020) 24318118
sales@sanner-in.com
www.sanner-group.com



AdCap®



AdCap® Desiccant Capsules

The AdCap® product family consists of drop-in desiccant capsules for the protection of moisture-sensitive drugs. This safe desiccant solution offers moisture protection on all sides: it adsorbs through a cardboard disc on the top as well as through a unique grid structure in the sidewall. The specially designed capsule ensures moisture adsorption wherever it is located in the container. What's more, the unique shape of AdCap® is designed to ensure maximum user protection.

Better moisture adsorbing performance

AdCap® stands for an innovative design that allows effective 360° moisture adsorption, regardless of the position in the container.

AdCap® does not lose its efficacy when positioned against the bottle sidewall or bottom.

Combined benefits of conventional capsules and canisters

AdCap® has a strong squeeze stability.

AdCap® features an innovative grid structure for improved moisture adsorption.

Non-dusting desiccants

AdCap® is non-dusting, because the capsule is filled with spherical silica gel or molecular sieve.

Patient safety

AdCap® reduces the risk of unintentional ingestion thanks to the unique grid structure.

Global standards

AdCap® is a cost-effective slot-in solution, developed according to global market needs and available all over the world.

AdCap® is compatible with standard dispensing equipment and filling lines.

Easy handling of packaging units

AdCap® comes in easy-to-handle quantities (packed in aluminium bags and carton), thus generating less weight and reducing residues per batch.

Readability

AdCap® offers an improved readability with a clear printing quality on the cardboard for safety indications.